

## Laboratory Accessories.



### TEST TUBE AND BURETTE BRUSHES


These brushes are ideal for cleaning cylinders, large test tubes, and bottles. Natural black bristle with splayed end provides a thorough cleaning action. The bristles are mounted in twisted galvanized wire. .

### TEST TUBE BRUSHES

D. Bristle, on galvanized wire. Head size 65x40mm Overall length 230mm

E. Tapered bristle with cotton tip, for tubes 19 to 25mm, bore. Overall length 230mm on galvanized wire.

F Bristle, for shaped end on galvanized wire, Head size 50 x 60mm. Overall length 230mm.




### BRUSH BEAKER

Nylon, head size 50 x 19 mm on galvanized wire, overall size 230 mm.



### BRUSH FLASK

Bristle head on wire handle.



### BRUSH BOTTLE

Hard nylon on iron wire. Size of brush 75x65 mm. Length of stem 375mm

### BRUSH BOTTLE

Nylon head on wire handle. & Bristle head on wire handle.

Head length x dia Overall For bottle

Length size

65 x 45mm 300 mm 125ml

75x50 mm 350mm 250ml

100 x 65mm 400 mm 500 ml

