



Works: 14, Pooja Vihar, Kardhan Road, Ambala Cantt
+91-171-3293762(O) 2699559(Fax), 9315100690(M)

Physics Instruments

1092-pH PAPER

Each vial contains 100 strips, 2 x 1/4; pH 1-14.
100 Strips



1093-SPATULAS

BOTH ENDS FLAT AND ONE END SPOON
STAINLESS STEEL SPATULA (Both ends flat)

B)Spatulas, with Wooden Handles
Flexible stainless steel blade with rounded end,
mounted in a polished wood handle. 3 blade, each



1094-TEST TUBE BRUSHES

D. Bristle, on galvanized wire. Head size 65x40mm Overall length 230mm
E. Tapered bristle with cotton tip, for tubes 19 to 25mm, bore. Overall length 230mm on galvanized wire.
F Bristle, for shaped end on galvanized wire, Head size 50 x 60mm. Overall length 230mm.



1095-BRUSH BEAKER

Nylon, head size 50 x 19 mm on galvanized wire,
overall size 230 mm.

