



Works: 14, Pooja Vihar, Kardhan Road, Ambala Cantt
+91-171-3293762(O) 2699559(Fax), 9315100690(M)

Physics Instruments



1112-CALORIMETER

WITH WOODEN BOX

Copper , 75 x 50 mm., height x diameter, resting on felt pad inside wooden polished box, which is provided with a brass nickel plated thermometer holder, with copper stirrer.



1113-CALORIMETER JOULE'S

For the determination of the specific heat capacity of liquid by the electrical method. The apparatus comprises a nickel-plated copper calorimeter 75 x 50 mm. fitted with a cork and a heating coil of thin constantan wire. The resistance of the coil is 6 ohms and the recommended working current is 0.5 A with a maximum of 1A. Supplied without Thermometer.



1114-PLOTTING COMPASS

With two clear glass faces 20mm diameter



1115-POCKET COMPASS

A small compass in a metal case, 45mm dia., marked in degree with N & S indications.

