



HIGH-QUALITY PHYSICS INSTRUMENTS

Bunsen Burner 183
Bunsen Burner, With Stopcock



Aneroid Barometer Demonstration 184
This is 100 mm dial graduated from 28" to 31" and 940 to 1040 millibars, demonstration type, with an index pointer, in transparent plastic case with pressure pump. It's only for educational demonstration purpose

Alcohol Burner, Aluminium 185
Alcohol Burner, Glass



Alcohol Burner Stand 186
Ideal for use with alcohol burners. Made of nickel plated steel wire gauge. Height 125mm.



Carbon Electrode 188

Burner, Methylated Spirit 187
Methylated Spirit burner is Metal, complete with wick, wick holder and cap.



Double Burette Clamp 189
Made of powder coated mild steel strip. The two clamps are provided with a spring action to hold the burettes. Fitted with a boss head for holding to a support rod.