



Works: 14, Pooja Vihar, Kardhan Road, Ambala Cantt
+91-171-3293762(O) 2699559(Fax), 9315100690(M)

Physics Instruments

1044-GLASS BLOCK

SEMI-CIRCULAR

Clear glass, all faces fully polished.
Dia x thickness




1045-GLASS BLOCK

RECTANGULAR

Made from moulded glass with polished surfaces. All angles accurately worked. For refraction experiment.

Size : 75x50x18mm, 100x60x25mm



1046-SET OF ACRYLIC BLOCK

Clear acrylic blocks 25 mm. thik. All faces fully polished comprise 7 block in storage tray as follows.
Rectangle 75 x 50 mm
Semi Circular 75 mm. dia.
Triangle 600 x 600 x 600 , 58 mm. side.
Triangle 900 x 450 x 450 , 75 „,hypotenuse.
Triangle 900 x 60 x 300 , 75 „,hypotenuse
Bi-Convex 75 mm. long, Covered faces 100mm.radius



1047-HOLLOW GLASS PRISMS

CUBES AND TRUFFS

From ordinary glass plates properly cemented with optical cement

Size 38 x 38mm Hollow Prism
Size 50x 50mm Hollow Prism
75x 50x 18mm Hollow Slabs
60 x 60 mm Hollow Truff
75 x 75mm Hollow Truff

