



Works: 14, Pooja Vihar, Kardhan Road, Ambala Cantt
+91-171-3293762(O) 2699559(Fax), 9315100690(M)


Physics Instruments



1133-INCLINED PLANE

WITH ANGLE MEASURER(metallic base)

A metal based board measuring 45 cm x 15 cm, hinged at one end supporting a scaled, planed board. The planed section can be inclined and fixed at any angle from 0-45 deg to the horizontal measured on the attached angle measurer. Without weights.



1134-STROBOSCOPE

A hardboard disc with twelve 50 x 3mm, slots and finger hold, pivoted on a ball bearing mount which also acts as a convenient handle. Diameter of disc 255mm.



1135-Prismatic Compass

Full Brass & Gun Metal Superior quality highly accurate and precisely manufactured fitted Superior fittings and complete with Tripod Stand Wooden or Aluminum

size 4", 4½", 5", 6".



1136-HOPE'S APPARATUS

To show the maximum density of water. Sheet iron cylinder mounted on a base and encircled midway along its length by a gallery 63 x 100 mm height x diameter. Fitted with tubulure, to carry thermometers and rubber stoppers.

